Nitrile Top

ref.no.	size
15390	6,5 - 7/S
15391	7,5 - 8/M
15392	8,5 - 9/L
15393	9,5 - 10/XL
15394	10,5 - 11/XXL

- DE Kategorie III - Chemische und biologische Risiken - FR Catégorie III - Risques chimiques et biologiques - EN Category III - Chemical and biological risks - IT Categoría III - Rischii chimici e biologici - NL Categorie III - Riesgos químicos y biológicos - ES Categoría III - Riscos químicos e biológicos - PT Categorie III - Chemiche en biologische risico's - SV Kategori III - Kemiska och biologiska risker - FI Luokka III - Kemialliset ja biologiset vaarat - DA Kategori III - Kemiske og biologiske risici - NO Kategori III - Kemiske og biologiske risikoer - PL Kategoria III - Zagrożenie chemiczne i biologiczne - HU III. Kategória - Kémiai és biológiai kockázatok - SK Kategória III - Chemicke a biologické rizíky - CS Kategorie III - Chemická a biologická rizika - SL Kategorija III - Kemična in biološka tveganja - HR Kategorija III - Kemikalije i biološke opasnosti - RU Kатегория III - Химические и биологические риски - RO Kategori III - Kimya-sal ve biyolojik riskler - TR Categoria III - Riscuri biologice și chimice

ISO 374-1
Type B

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Instructions for use

The glove types named here meet with the requirements of regulation (EU) 2016/425, EN ISO 374-1:2016 and EN 420:2003+A1:2009. As the gloves are adjusted to special purposes, the lengths may deviate from the requirements of EN 420:2003+A1:2009.

Results chemical tests:

n-heptano (I):	Nivel 3
40% Sodium Hydroxide (K):	Nivel 6
37% Formaldehyde (T):	Nivel 6

Level	1	2	3	4	5	6
Breakthrough time (min)	> 10	> 30	> 60	> 120	> 240	> 480

EN 374-4:2013 - Resistant to degradation against chemicals:

n-heptano (I):	38,3 % (+/- 4,8 %)
40% Sodium Hydroxide (K):	-6,8 % (+/- 4,4 %)
37% Formaldehyde (T):	23,7 % (+/- 9,9 %)

Warning / risk assessment: a) Gloves to protect against mechanical action whose effects are superficial, substances and mixtures which are hazardous to health, and harmful biological agents. Important: The gloves are recommended for use in situations where only low chemical protection is identified as needed. b) While selecting an equipment, user should perform risk analysis based on the intended use and determine the suitability of the product's test standards and protection levels obtained. c) Information provide does not reflect the actual duration of protection in the workplace due to other factors influencing the performance, such as temperature, abrasion, and degradation, and the differentiation between mixtures and pure chemicals. d) Information regarding protection refers to the working surface, i.e. 'the palm' of the glove, which has been submitted to testing. e) Gloves should be thoroughly inspected for damages (specify for nicks and holes) before use. f) The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. g) It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. h) When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. i) The maximum wear time depends on the activity being carried out and the person. j) EN374-4:2013 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. k) The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. l) This product contains Nitrile Butadiene Rubber and compounding chemicals which may cause allergic reaction in some individuals. m) The glove offers no protection against perforation with sharp objects, e.g. injection needles.

Storage: To be stored in cool dry place away from sunlight. Storage procedures are the main factor in determining glove shelf life. Gloves should be kept in their packaging protected from sunlight, artificial light, humidity and stored at temperatures between 5 °C - 30 °C. Real time stability testing for expiration date on products after accelerated ageing process are declared as 5 years.

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